



A-LINETM

Designed specifically for aspiration of blood sample from arterial lines.



BD LUER-LOCKTM SYRINGE

- Recommended for a secured connection of Syringe to any device (Needle, tip cap etc.)
- Prevents unexpected disconnection during Sampling and Handling.



BD HEMOGARD™ - TYPE TIP CAP

- Minimum Contact with Blood
- Screw Lock ensures tight seal even in transportation
- Large size cap ensures easier handling and designed to prevent blood contact with contaminated syringe tip



QUALITY SYRINGE WITH SPRAY DRIED BALANCED LITHIUM HEPARIN

Eliminating Preanalytical errors for the best patient outcome

CERTIFICATIONS



Registered with Drug Regulatory Authority of Pakistan (DRAP)



U.S FDA Approved: US Food & Drug Regulatory Authority Requires medical devices to demonstrate safety & efficacy prior to marketing. BD ABG Syringes are made on FDA approved plant which ensures both parameters



ISO 11137 Certified: Which ensures Gamma Radiation sterilization - highest quality sterilization method

AVAILABLE SIZE

Ref#	Description	Balanced Lithium Heparin	Tip	Tip Cap	Case Qty
364378	3ml Aspiration Syringe	80 I.U.	Luer Lock	BD Hemogard	100

FEATURES & BENEFITS

FEATURE	ADVANTAGE	VALUE	
Pre-Heparinized Syringes	Ensures correct and consistent concentration of Heparin according to IFCC guidelines Minimal chances of a clotted sample	Perfect sample for Blood Gas Reduces repeat samples requests	
Lithium Heparin	No interference in Sodium estimations	Lesser dependence on clinical judgments Reduces risk of falsely elevated Sodium value	
Calcium Balanced Heparin	Binding site of electrolytes in heparin saturated	Accuracy in measurements of electrolytes like Potassium / Calcium since they are relatively low in concentration	
Spray Dried Heparin	No dilution of sample	Accuracy in small volume samples Hb estimation is more accurate	
Closure Cap & High-Density Polypropylene Material	Helps maintain Anaerobic state of sample No need to bend the needle High Density Polypropylene ensures no gaseous exchange.	No contact with air while transportation of sample Safety for user	

Imported & Marketed by:
ASTO Life Sciences Pvt. Ltd

Address: Plaza # 1, Block Orchard 1, Paragon City, Barki Road , Lahore
Email: info@astolife.com
website: www.astolife.com

